



DENTAL DRILL REPAIR



CuttingEdge

Surgical Instruments Ltd

Cutting Edge Surgical Instruments Ltd offer a full and comprehensive repair service for all makes and models of dental handpieces and scaler equipment.


All handpieces are repaired to OEM specifications. The basic service comprises of a full strip-down, with all parts being ultrasonically cleaned and checked for wear, all bearings and seals replaced, collets reset, reassembled, and tested.

Our fully trained technicians have a wealth of experience with handpieces from all the leading brands.

Prices

Tier 1 Handpiece repair	£TBA
Tier 2 Handpiece repair	£TBA
Irrigation nozzles only	£TBA

All prices subject to change.



CuttingEdge

Surgical Instruments Ltd

T · 01992 558 222

E · info@cuttingedgesurgical.co.uk

www.cuttingedgesurgical.co.uk

4B Great Northern Works,
Hartham Lane,
Hertford, SG14 1QN