

AIR TOOL REPAIR SERVICE



CuttingEdge

Surgical Instruments Ltd

AIR TOOL REPAIR SERVICE

Cutting Edge Surgical Instruments Ltd offer a comprehensive repair service for all makes and types of Air Tools, Hoses and associated Attachments.

All repairs are undertaken by fully trained and experienced technicians.

All handpieces are repaired to the OEM specification.

Our service comprises of a full strip down.

All parts are ultrasonically cleaned and checked for wear.

All bearings, seals, and air veins are replaced.

All moving parts are re-lubricated.

Handpieces are then re-assembled and undergo rigorous testing.

We also reset any collets and chucks to original specifications and replace inner hoses if required.

We repair and service items from all major manufacturers including Howmedica, Microaire, Stryker, Synthes, 3M, Zimmer/Hall, and Desoutter.

We collect your repairs and deliver them back to you free of charge.

All repairs are returned with a technical report detailing work carried out and any parts used.

All repairs are fully covered by a one year guarantee against faulty workmanship and parts.

For more information and to arrange a collection please contact us.



CuttingEdge

Surgical Instruments Ltd

... a cut above the rest!

tel: 01992 558 222

email: info@cuttingedgesurgical.co.uk

website: www.cuttingedgesurgical.co.uk

Unit 1 Harforde Court, John Tate Road

Hertford, Hertfordshire SG13 7NW