


BATTERY POWER TOOL REPAIR



CuttingEdge

Surgical Instruments Ltd

BATTERY POWER TOOL REPAIR

Cutting Edge Surgical Instruments Ltd offer a comprehensive repair service for most makes and types of Battery Power Tools.

All repairs are undertaken by fully trained and experienced technicians.

All Power Tools are repaired to the OEM specification.

Handpieces will be performance tested, fully stripped down, comprehensively interior cleaned, bearings and o-rings replaced, controller checked and serviced, motor tested and serviced, lubricated, then finally rebuilt and tested.

We collect your repairs and deliver them back to you free of charge.

All repairs are returned with a technical report detailing work carried out and any parts used.

All repairs are fully covered by a one year guarantee against faulty workmanship and parts.

We can continue to repair and service many older handpieces even when the OEM claims obsolete.

High Quality Repairs, Fast Turnaround Times, and Competitive Prices.



CuttingEdge

Surgical Instruments Ltd

... a cut above the rest!



tel: 01992 558 222

email: info@cuttingedgesurgical.co.uk

website: www.cuttingedgesurgical.co.uk

Unit 1 Harforde Court, John Tate Road

Hertford, Hertfordshire SG13 7NW