

DENTAL DRILL REPAIR SERVICE



CuttingEdge

Surgical Instruments Ltd

DENTAL DRILL REPAIR SERVICE

Cutting Edge Surgical Instruments Ltd offer a comprehensive repair service for all makes and types of Dental Handpieces.

All repairs are undertaken by fully trained and experienced technicians.


All handpieces are repaired to the OEM specification.

Handpieces will be fully stripped down, all parts ultrasonically cleaned and checked for wear, all bearings and seals are replaced, collets reset, reassembled and tested.

We collect your repairs and deliver them back to you free of charge.

All repairs are returned with a technical report detailing work carried out and any parts used.

All repairs are fully covered by a one year guarantee against faulty workmanship and parts.



All makes and models of handpieces can be serviced and repaired.

High Quality Repairs, Fast Turnaround Times, and Competitive Prices.



CuttingEdge

Surgical Instruments Ltd

... a cut above the rest!

tel: 01992 558 222

email: info@cuttingedgesurgical.co.uk

website: www.cuttingedgesurgical.co.uk

Unit 1 Harforde Court, John Tate Road

Hertford, Hertfordshire SG13 7NW