

# ELECTROSURGICAL CABLES

**CuttingEdge**

Surgical Instruments Ltd

## BIPOLAR, MONOPOLAR AND LAPAROSCOPIC CABLES

Machine Connector

Product Description

Instrument Connector



Bipolar Lead with 2 x 4mm Pins and Standard Euro Instrument Block



Monopolar Lead with 4mm Machine End and 4.8mm Socket Instrument End



Monopolar Lead with 8mm Machine End and 4.8mm Socket Instrument End



Laparoscopic Lead with 4mm Machine End and 4mm Socket Instrument End



Laparoscopic Lead with 8mm Machine End and 4mm Socket Instrument End

# ELECTROSURGICAL CABLES

Cutting Edge Surgical Instruments Ltd supply a range of Monopolar and Bipolar Cables for regular and Laparoscopic Surgical Instruments and Devices.

Our Reusable Monopolar Diathermy and Bipolar Cables are high quality electrosurgical cables made to supply HF current from the electrosurgical generator to Monopolar and Bipolar Forceps and Electrosurgical Accessories.

These are high quality and reusable cables made using superior materials ensuring optimum performance and issue free use.

Our cables are available with all popular generator and instrument connectors including 4mm sockets for Laparoscopic Instruments.

We keep very large stocks of all cables for fast deliveries.

For more information and to place an order please contact us.



**CuttingEdge**

Surgical Instruments Ltd

*... a cut above the rest!*



tel: 01992 558 222

email: [info@cuttingedgesurgical.co.uk](mailto:info@cuttingedgesurgical.co.uk)

website: [www.cuttingedgesurgical.co.uk](http://www.cuttingedgesurgical.co.uk)

Unit 1 Harforde Court, John Tate Road  
Hertford, Hertfordshire SG13 7NW