

PHACO HANDPIECE REPAIR



CuttingEdge
Surgical Instruments Ltd

PHACO HANDPIECE REPAIR

Cutting Edge Surgical Instruments Ltd offer a comprehensive repair service for most makes and types of Phaco Handpieces.

All Handpieces will be fully inspected, dismantled and fully stripped, we replace the piezoelectric crystals, connector, cord, and strain relief, new o-rings, replacement wiring, new chipset and shims.

Then we rebuild and test the handpieces to ensure it works again like a brand new handpiece.

We can save you up to 70% of the cost of a repair from the OEM.

We collect your repairs and deliver them back to you free of charge.

All repairs are returned with a technical report detailing work carried out and all parts used.

All repairs are fully covered by a one year guarantee against faulty workmanship and parts.

Please contact us for more information.



CuttingEdge

Surgical Instruments Ltd

... a cut above the rest!

tel: 01992 558 222

email: info@cuttingedgesurgical.co.uk

website: www.cuttingedgesurgical.co.uk

Unit 1 Harforde Court, John Tate Road

Hertford, Hertfordshire SG13 7NW