

SET RECONDITIONING



Bringing Sets Back to
an 'As New' Condition

CuttingEdge

Surgical Instruments Ltd

SET RECONDITIONING

Cutting Edge Surgical Instruments Ltd offer a fast and reliable Set Reconditioning service. Our Set Reconditioning service ensures your instrument sets are returned to an 'As New' condition. Every item on the set will be inspected and then fully reconditioned, repaired, serviced or replaced. All Diathermy and Bipolar Forceps will be stripped, recoated, and rebuilt and leads replaced. All Needleholders will have new jaws fitted and gold handles replated. All Scissors will cut correctly again. All Forceps will grip again. All Retractors will retract again. Every repaired and serviced item will be etched with the date of repair for traceability purposes. All repairs are fully covered by a one year guarantee against faulty workmanship and parts. All sets are collected and returned to you free of charge. Part of the Set Reconditioning process is the microscopic inspection of every instrument on the set. Inspection and preventative maintenance will ensure that the minor issues we find and repair do not become major issues with the potential of rendering the instrument beyond repair, and requiring a costly replacement, or even causing incidents in surgery.


Please contact us for more information.



CuttingEdge

Surgical Instruments Ltd

... a cut above the rest!



tel: 01992 558 222

email: info@cuttingedgesurgical.co.uk

website: www.cuttingedgesurgical.co.uk

Unit 1 Harforde Court, John Tate Road

Hertford, Hertfordshire SG14 1QN